

HelixAmp™ FastLAMP Kit (Ver. 2.0)

Kit Contents

HelixAmp™ FastLAMP Kit (Ver. 2.0)		
Cat. No.	FLMP2-100 (100rxns)	FLMP2-500 (500rxns)
FastLAMP Enzyme V2	0.1ml	0.1ml x 5ea
5x FastLAMP Buffer V2 (Mg-free)	0.5ml	0.5ml x 5ea
100mM MgSO ₄	0.25ml	0.25ml x 5ea
D-Solution	1ml	1ml x 5ea
Instructions for Use	1ea	1ea


Description

HelixAmp™ FastLAMP Kit (Ver. 2.0) provides simple and fast (within 20 minutes) target DNA amplification using loop-mediated isothermal Amplification (LAMP). This kit consists of 5x FastLAMP buffer V2 (Mg-free), FastLAMP Enzyme V2 and a D-Solution. The 5x FastLAMP Buffer V2 (Mg-free), optimized for fast amplification, contains buffering reagents, dNTPs, and salts. The FastLAMP Enzyme V2 is composed of a newly engineered *Bst* DNA Polymerase that provides improved amplification reaction properties. The novel *Bst* DNA polymerase enhances the DNA polymerization speed and allows the fast isothermal amplification reaction.

Application

Isothermal amplification (LAMP) of DNA target

Storage



Store below -20°C

Shelf life



12 months

Quality Control

By NanoHelix's ISO 13485-certified Quality Management System, each lot of **HelixAmp™ FastLAMP Kit (Ver. 2.0)** was tested against predetermined specifications to ensure consistent product quality.

Protocol

1. Reaction Mixture

- LAMP products can be analyzed by examining the end-point product or real-time assay.
- Prepare the reaction mix according to the following table for the selected analysis method.

Components	For end-point	For real-time (intercalating dye or probe)
Template ¹⁾	1 ~ 5 μ l	1 ~ 5 μ l
5x FastLAMP Buffer V2 (Mg-free)	5 μ l	5 μ l
FastLAMP Enzyme V2	1 μ l	1 μ l
100mM MgSO ₄	2.25 μ l	2.25 μ l
10x LAMP Primer mix ²⁾	2.5 μ l	2.5 μ l
Fluorescent dye ³⁾ or probe	-	x μ l
RNase-free Water	Adjust to final 25 μ l	

¹⁾ **Template preparation:** Purified DNA sample using a commercial DNA preparation kits or a manual method can be applied directly to this assay. For better amplification and detection sensitivity, we recommend to use the **D-Solution** provided in this kit. **D-Solution** helps to denature the template DNA and induces efficient primer binding to its target sequence. The denatured template can be prepared by adding 1/10 volume of the D-Solution to the DNA sample.

Ex) DNA Sample 50 μ l + D-Solution 5 μ l

²⁾ For simplicity in setting up reactions, we recommend making stocks of the LAMP primers at a usable concentration. For example, we suggest a following **10x LAMP Primer Mix** containing all six LAMP primers. If there is low-efficiency or non-specific amplification, modify the primer concentration or design a new set of primers for the target sequence.

10x LAMP Primer Mix	
LAMP primers	Primer concentration.
FIP	16 μ M
BIP	16 μ M
F3	2 μ M
B3	2 μ M
LF	8 μ M
LB	8 μ M

- 3) Recommend using final 0.1~0.2x SYBR Green I or 0.1~0.3x EvaGreen dye (not supplied in this kit).

2. Reaction Condition

For end-point assay: Incubate at 65°C for 30 minutes. Time can be extended as necessary for very low copy targets, challenging sample types, etc. Analysis the reaction product by a gel-electrophoresis or other detecting tools including colorimetric and fluorescence detection, turbidity observation, lateral flow devices, etc.

For real-time assay (intercalating dye or probe): Use a real-time PCR machine or an isothermal amplification instrument to run the assay. Set the instrument to a constant incubation temperature at 65°C. Measure the fluorescence intensity at every 1 min for 30 minutes. The reaction time can be increased as necessary for very low copy targets, challenging sample types, etc.

Products

Cat. No.	Products	Size
FLMP2-100	HelixAmp™ FastLAMP Kit (Ver. 2.0)	100rxns
FLMP2-500	HelixAmp™ FastLAMP Kit (Ver. 2.0)	500rxns